# Thopaz™

## DIGITAL CHEST DRAINAGE SYSTEM

- EN Instructions for use
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- DA Brugsanvisning
- ES Instrucciones de uso
- PT Instruções de utilização
- PL Instrukcja obsługi


# Thopaz University

At the web address below you will find short how-to-use videos for Thopaz.

[www.medela.com](http://www.medela.com)

Thopaz University is designed to support you. Each video lasts no longer than 60 seconds and includes text bullet points and voice commentary.

Thopaz University can be found at [www.medela.com](http://www.medela.com) via the “University” tab. It is available on all formats including desktop, smartphone and tablet. We recommend saving Thopaz University as a shortcut on desktop, smartphone and tablets for even quicker access.

# Table of Contents

- Congratulations ..... 6**
- 1 Warnings and safety instructions ..... 7**
  - Safety instructions ..... 9
- 2 Power supply and battery operation ..... 10**
  - Backup battery ..... 10
- 3 Charge battery ..... 11**
  - Battery low ..... 11
  - Battery empty ..... 11
- 4 Description ..... 12**
  - Introduction ..... 12
  - Intended use/indications ..... 12
  - Contraindications ..... 12
  - Intended user ..... 12
  - Intended patient population ..... 12
  - Important note ..... 12
- 5 Overview ..... 13**
  - Definition of vacuum ..... 13
  - Main elements of Thopaz ..... 13
  - Operating elements ..... 14
  - Activating buttons ..... 14
  - Tubings ..... 15
  - Canisters ..... 17
- 6 Additional information ..... 18**
  - Regulated pressure ..... 18
  - Safety pressure limit ..... 18
  - Tube flushing ..... 18
  - Safety chamber in canister ..... 18
  - Overflow protection/bacteria filter ..... 18
  - Leakage warning (does not correspond to «Leak in system» alarm) ..... 18
- 7 Installation ..... 19**
  - Check initial delivery ..... 19
  - Initial startup ..... 20
- 8 Preparation for use/operating instructions ..... 22**
  - Checks before use ..... 22
  - Connect tubing ..... 22
  - Snap in canister ..... 23
  - Using 2L canister ..... 23
  - Switch Thopaz on ..... 24

Carry out functional check .....	25
Check therapy progress.....	26
Rounding Function .....	26
Graph mode.....	27
Zoom mode.....	28
Catheter check.....	28
Switching Thopaz off.....	29
Transferring data to PC with ThopEasy+ .....	29
<b>9 Replace canister .....</b>	<b>30</b>
Canister with solidifier: activation of the solidification process.....	31
<b>10 Changing pressure.....</b>	<b>32</b>
Change pressure.....	32
Gravity mode.....	32
<b>11 Changing settings .....</b>	<b>33</b>
Change settings .....	33
Factory settings.....	34
Setting a changed parameter as standard.....	34
<b>12 Taking a drainage sample .....</b>	<b>35</b>
<b>13 Troubleshooting and alarms .....</b>	<b>36</b>
Pump not running.....	36
Description .....	36
Canister full.....	38
<b>14 Warranty, maintenance and checks.....</b>	<b>39</b>
Warranty.....	39
Maintenance.....	39
Routine Check.....	39
Instructions for FPT test .....	40
Safety Check.....	41
<b>15 General cleaning guidelines .....</b>	<b>42</b>
General notes .....	42
Disposables.....	42
<b>16 Accessories .....</b>	<b>44</b>
<b>17 Signs and symbols.....</b>	<b>45</b>
<b>18 Technical specifications .....</b>	<b>47</b>
<b>19 Disposal .....</b>	<b>48</b>
<b>20 Technical documentation .....</b>	<b>473</b>
<b>Appendix A.....</b>	<b>478</b>
<b>Appendix B.....</b>	<b>480</b>

# Congratulations

With Thopaz, you have acquired an innovative digital chest drainage system that sets new standards in therapy and chest drainage management.

Thopaz has an electronic measuring and monitoring system with optical and acoustic status indications. The device is a dry system, which means no fluids are necessary for operation. Important information regarding the course of therapy is indicated digitally and graphically on the display. Upon completion of therapy, it can be transferred to a PC.

The handy size and extreme quietness of operation are pleasant for the patient as well as those in the vicinity.

# 1 Warnings and safety instructions

## **WARNINGS**

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

## **CAUTIONS**

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

## **Safety related tip**

Indicating useful information about the safe use of the device.

Thopaz is approved exclusively for the use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when Thopaz is used in combination with the original Medela accessories (canisters, tubings, carrying strap, mains adapter, docking station – see chapter 16).

Thopaz is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfil the requirements as outlined in the IEC 60601-1-2 standard. Untested HF sources, radio networks or the like can impair the function of the device and should not be operated in the vicinity of Thopaz.

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1-1 or clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local representative or the technical service department.

**Please read and observe these warning and safety instructions before operation. These instructions for use must be kept with the device for later reference.**

Please note that these instructions for use are a general guide for the use of the product. Medical matters must be addressed by a physician.

Medela considers itself only responsible for the effect on BASIC SAFETY, reliability and performance of Thopaz if it is used in accordance with the instructions for use.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Subject to change.

## WARNINGS

- Consult this manual prior to use.
- Data transfer via USB is not allowed during therapy.
- Before cleaning the device, pull the plug out of the fixed mains socket.
- Do not use Thopaz in MRT (Magnetic Resonance Tomography).
- Do not dry Thopaz with microwaves.
- The device is not suitable for use while bathing, showering or in a hazardous explosive environment.
- The pressure range to be set must be determined by a physician in accordance with the age and weight of the patient.
- The secretions gathered must be monitored more regularly or frequently in accordance with the age and weight of the patient.
- To ensure Thopaz can reduce the pressure, it is essential to wait 30 seconds between taking the sample from the sample port and unclamping the patient catheter (drain).
- Do not use Thopaz if the drainage therapy indicates a pressure greater than Thopaz max. pressure range of -10 kPa.
- Do not use Thopaz if the drainage therapy indicates a flow rate greater than Thopaz max. flow capacity of 5 l/min.
- Do not use Thopaz if the drainage therapy indicates no pressure should be applied to the patient.
- If a persistent air leak is shown, check the system is assembled correctly before taking further corrective actions. Ensure the system is air tight by clamping the catheter and observing the air leak decrease to zero.
- The catheter/connector interface is a location where clotting may occur. We recommend regular monitoring of this interface and appropriate removal procedure in case of occurrence.
- Do not connect bilateral thoracic drains to one Thopaz unit. In such cases, the use of two Thopaz units is recommended.
- The use of Thopaz for any other indication than intended is neither desired nor allowed.
- The Thopaz system is not designed for retransfusion.
- This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire.
- For pediatric patients adapt pressure settings according to hospital guidelines.
- Do not touch the patient and the docking contacts simultaneously.

## CAUTIONS

- No modification of all articles related to the Thopaz system is allowed.
- Before connecting Thopaz to the power supply, please verify that the supply voltage corresponds to that given on the device specification plate.
- Incorrect use can cause pain and injury to the patient.
- Do not use sterile accessories when the sterile packaging is damaged.
- For patients for whom a breakdown of the device can lead to a critical situation a replacement device must always be available.

## Safety related tip

- A visual inspection of the housing and power cord must be carried out before each use.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.

## Safety instructions

- Thopaz is a medical device that requires special safety measures in regard to EMC. It must be installed and put into operation in accordance with the attached EMC information in chapter 22.
- Portable and mobile RF communications equipment can affect medical devices.
- In each of the following cases, Thopaz must not be used and it must be repaired by Customer Services:
  - if the power cord or the plug are damaged.
  - if the device is not functioning according to manufacturer.
  - if the device is damaged.
  - if the device shows safety defects.
- Keep the power supply cord away from hot surfaces.
- The mains plug and the on/off switch must not come into contact with moisture. Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord!
- Thopaz must stand upright during use.
- Separation from the mains is only assured through the disconnection of the mains adapter and fixed socket connection.
- Never use the device at high room temperatures, while bathing or showering, if you are very tired or in an environment where there is a risk of explosion.
- Never place Thopaz in water or other liquids.
- When using single use, sterile products, please note that they are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or

- biological characteristics. Reuse could cause cross contamination.
- Contact your local Medela customer service representative for assistance with product operations.
  - Wear gloves for all operations.
  - The patient should be regularly monitored according to internal hospital guidelines.
  - A replacement device must always be available for patients for whom a breakdown of the device can lead to a critical situation.

These instructions for use must be kept for later reference.

## 2 Power supply and battery operation

### WARNINGS

This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire.

Thopaz can be operated from the mains power supply [REF](#) 077.0148, from the docking station [REF](#) 079.0037 or with the integrated lithium-ion battery. Before using the docking station, read first the docking station instructions [REF](#) 200.1554. The battery is charged on mains operation or on docking station. The battery duration is dependent on the run-time of Thopaz. This is influenced by the extent of parenchymal leakage and the set pressure. Thopaz does not run continuously but only switches on when the actual and nominal negative pressure differ. During continuous operation, Medela guarantees a minimum run-time of 4 hours in battery operation after the device has been fully charged. In practical operation the real battery run time results in > 10 hours.

If Thopaz is not used very often, the battery must be charged approximately every 6 months to ensure optimum functioning.


### Safety related tip

#### Backup battery

If an internal defect occurs (broken cable, battery defective), an acoustic warning signal is sounded for at least 3 minutes. Power is supplied by the backup battery. Replace Thopaz immediately.

# 3 Charge battery

Plug the Thopaz Power supply into a wall outlet and open the small protection cover to access the DC socket. Plug the coaxial power supply into the DC socket. Moving bars on the battery display will indicate charging is occurring [ ].

## Battery low

Approximately 30 minutes before the battery is fully empty an acoustic signal sounds and the battery symbol [ ] starts to blink. The acoustic signal can be muted but the battery symbol continues to blink until the battery is empty. The set pressure is maintained but the battery should be recharged as soon as possible.

to acknowledge an acoustic alarm (see chapter 13)

Press simultaneously, then follow the instructions on the display

## Battery empty

The battery symbol will blink with an acoustic sound for 10 minutes prior to the battery fully discharging, unless the pump is switched off earlier. The acoustic sound cannot be muted during the 10 minutes. If Thopaz is not connected to a power source the system will switch off automatically after 10 minutes. The negative pressure is not maintained after the system turns off.

1. Battery empty (sign is blinking)

3. Battery is being charged (moving bars)

2. Battery almost empty

4. Battery is fully charged

# 4 Description

## **Introduction**

Thopaz is a high-quality digital chest drainage system. The compact system provides regulated negative pressure close to the patient's chest. It monitors the air leak and only applies the suction required to maintain the negative pressure prescribed and set by the managing physician. The digital display provides objective data in real time as well as a 24-hour historical graph of the air leak, which allows easy tracking of the therapy progress. The system combines easy handling and cleaning with safety features to ensure optimal operation.

## **Intended use/indications**

Thopaz is intended to be used for aspiration and removal of surgical fluids, tissue, gases, bodily fluids or infectious materials. Thopaz is indicated for all situations where chest drains are applied - especially for thoracic drainage in the pleural and mediastinal cavity in situations such as pneumothorax, after cardiac or thoracic surgery (post-operative), thorax injury, pleural effusion, pleural empyema or other related conditions. Thopaz is intended for use on patients in appropriate care settings.

## **Contraindications**

There are no known contraindications for the Thopaz digital chest drainage system.

## **Intended user**

Thopaz should only be operated by trained users. These persons must not be hard of hearing or deaf and must have adequate visual faculty.

## **Intended patient population**

Thopaz is intended to be used on patients only exhibiting conditions as described in the indications for use.

## **Important note**



Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.

# 5 Overview


## Definition of vacuum

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in Kilopascal (kPa). In this document, the indication of -10kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure (according to EN ISO 10079:1999).

## Main elements of Thopaz

## Operating elements

 **Press simultaneously**   
**to change pressure and settings**  
(see chapter 10, 11)

## Activating buttons

Press the buttons briefly to activate.

**Exception:** The «Standby» button must be pressed for longer than 3 seconds to prevent an unintended interruption of the therapy.



## Tubings

**STERILE EO**

Normal/small/large/single/double connector

Material: PVC (medical grade) sterile, twofold packaging

Length: 1.5 m

079.0025/26: Patient tube inner diameter: 5.5 mm  
 Connector outside diameter 1: 9.5 mm  
 Connector outside diameter 2: 16 mm

079.0021/22: Patient tube inner diameter: 5.5 mm  
 Connector outside diameter 1: 6.4 mm  
 Connector outside diameter 2: 16 mm

079.0023/24: Patient tube inner diameter: 5.5 mm  
 Connector outside diameter 1: 4.0 mm  
 Connector outside diameter 2: 11.0 mm

Thopaz tubing single,  
large connector

**REF** 079.0025

Thopaz tubing double,  
large connector

**REF** 079.0026

Thopaz tubing single

**REF** 079.0021

Thopaz tubing double

**REF** 079.0022

Thopaz tubing single,  
small connector

**REF** 079.0023

Thopaz tubing double,  
small connector

**REF** 079.0024

### **Safety related tip**

Medela recommends to use the biggest possible connector size compatible for the catheter.

**Canisters**

**STERILE EO**

With/without solidifier 0.3l/0.8l/2l


Material: Polypropylene, sterile

Reading accuracy: +/- 2.5% (in upright position)

Solidifier: 0.3l = 12g/0.8l = 28g/2l = 2x30g

Simethicone capsule: prevents foaming

(0.3 l = 720mg/0.8l = 1440mg/2l = 1440mg Simethicone)

Thopaz canister 0.8l, antifoaming **REF** 079.0016

Thopaz canister 0.8l with solidifier, antifoaming

**REF** 079.0017

Thopaz canister 0.3l, antifoaming

**REF** 079.0011

Thopaz canister 0.3l with solidifier, antifoaming

**REF** 079.0012

Thopaz canister 2l, antifoaming

**REF** 079.0018

Thopaz canister 2l with solidifier, antifoaming

**REF** 079.0019

**👉 Safety related tip**

Use sample port of tubing to drain a sample. Therefore refer to chapter 12 – Taking a drainage sample.

**👉** The simethicone capsule requires 60 min of submersion for full effect.

# 6 Additional information

## **Regulated pressure**

The pressure is regularly checked and adjusted by Thopaz. The set pressure corresponds to the pressure of the entire system including the patient.

## **Safety pressure limit**

The adjustable pressure range for Thopaz is -0,1 to -10kPa (-100 mbar, -75 mmHg). If pressure value is greater than -7 kPa (-70 mbar, -53 mmHg) are set, the following warning appears in the display: too high pressure can cause pain and injuries to the patient. This warning must be acknowledged with «OK» before the pressure can be increased further.

## **Tube flushing**

Flushing of the tubing occurs approximately every 5 minutes or when Thopaz detects a siphon and prevents clogging of the patient tubing.

## **Safety chamber in canister**

For appropriate, safe operation Thopaz must be in an upright position. If Thopaz tips, the special construction of the safety chamber in the upper region of the canister protects the hydrophilic filter from premature clogging. The pressure is maintained. Important: immediately set Thopaz upright again. This feature is disabled if the maximum filling capacity of the canister is exceeded.

## **Overflow protection/bacteria filter**

A hydrophilic overflow protection/bacteria filter in the canister as well as in the tube sets (in the connector to Thopaz) protects Thopaz from the entry of fluids and protects the environment from contamination.

## **Leakage warning (does not correspond to «Leak in system» alarm)**

A significant and sustained increase in the flow rate indicates that there is a leakage in the system or an irregularity in the course of therapy. This function can be switched on as required - see chapter 11.

# 7 Installation


## 7.1 Check initial delivery

Check the delivery package of Thopaz for completeness and general condition.

Thopaz™

**REF** 079.0000 with Instructions for use 200.0685  
079.0002 with Instructions for use 200.0686  
079.0003 with Instructions for use 200.2004 / 101034234

Thopaz wall mains adapter coaxial DC plug

**REF** 077.0148

Thopaz CD

**REF** 079.0035

Instructions for use

**REF** 200.0685 (for 079.0000)  
200.0686 (for 079.0002)  
200.2004 / 101034234  
(for 079.0003)

## 7.2 Initial startup

### WARNINGS

For use only after instruction by medically trained and qualified persons who have been adequately trained in the use of Thopaz.

### CAUTIONS


Before the first application in accordance with the intended use, an inspection according to IEC 62353 (see Appendix B) is recommended.

## 7.2a Checks before use

- Check Thopaz system before use for damage of the power cord or plug, obvious device damage or safety defects and proper functioning of the device.
- Check for completeness and general condition of Thopaz delivery package.

## 7.2b Initial startup procedure


1.

1.1 Select plug.

1.2 Attach plug to the mains adapter [REF](#) 077.0148.

2.

2.1 Connect Thopaz to mains by mains adapter [REF](#) 077.0148 or by docking station [REF](#) 079.0037.


**Important:** Do not yet connect the patient to the system.

3.


3.1 Charge for approximately 3 hours or until the [  ] symbol appears in the display.

3.2 Press [  ] to switch Thopaz on. A selftest starts.

4.  4.1 Select language.  
4.2 Confirm with «OK».  
4.3 Follow instructions on display.

5.  5.1 Press [  ] to switch Thopaz off.

6.  6.1 Disconnect Thopaz from the fixed mains socket by pulling on the plug housing. Do not pull on the cable or the bend protection.

7. Thopaz is ready for preparation for use.

# 8 Preparation for use/operating instructions

## WARNINGS

For use only after instruction by medically trained and qualified persons who have been adequately trained in the use of Thopaz. Wear gloves for all operations.  
«Pressure» generally implies «negative pressure».

## CAUTIONS

- Sterile accessories must be checked on the integrity of the packaging before use.
- Non sterile and reusable accessories must be clean and disinfected according to the Medela cleaning guide 200.2391
- Do not switch Thopaz on before charging the battery first.

## 8.1 Checks before use

- Check Thopaz system before use for damage of the power cord or plug, device damage or safety defects and proper functioning of the device.
- Check that the seal is correctly placed in the vacuum port of Thopaz. If you do not know where the vacuum port is placed please refer to Chapter 5 – Overview.
- Check for completeness and general condition of the Thopaz delivery package.
- Make sure that the rechargeable battery is recharged in case Thopaz is operated in battery mode.
- Check all accessories prior to use:
  - canister for cracks, brittle and flawed spots.
  - tubing for cracks, brittle areas and that connectors are firmly attached. Replace if necessary.


## 8.2 Connect tubing

1.

1.1 Open external packaging.

2.


2.1 Insert the smaller of the two connectors (connection to pump) horizontally in the direction of the arrow.

## 8.3 Snap in canister

1. 
  - 1.1 Unpack the canister.
  - 1.2 Position the openings at the top and attach the bottom part of the canister to Thopaz.
2. 
  - 2.1 Push the canister into Thopaz until you hear a click.


### Using 2L canister

1. 
  - 1.1 Unpack the canister.
  - 1.2 Position the openings at the top and attach the bottom part of the canister to Thopaz.
  - 1.3 Push the canister into Thopaz until you hear a click.

### Safety related tip

Never place Thopaz with an attached 2l canister on the optional docking station without attaching the optional Adapter docking station [REF 079.0038](#) for Thopaz Canister 2L first.


### CAUTIONS

The Thopaz 2l canisters are not intended for portable use (not by hand and not using the carrying strap).

## 8.4 Switch Thopaz on

**Important:** Do not yet connect the patient to the system.

1.
  - 1.1 Press [  ] to switch on Thopaz. A selftest starts.
  - 1.2 Check acoustic signal (beep).

If the selftest is not successful, see the troubleshooting instructions on the display and refer to chapter 13.

2.
  - 2.1 Confirm whether a new patient is connected or not. The therapy number is important for data transfer to the PC.

**yes** Thopaz will issue a new therapy number (recommended for a new patient).

**no** The therapy number remains unchanged (recommended for the continued treatment of the same patient). The graphic presentation is retained.

## 8.5 Carry out functional check

1. 
  - 1.1 Seal off the conical connector of the patient tubing with your thumb (through the packaging). The tubing must not be kinked otherwise the measuring tube will be closed.


2. 
  - 2.1 Switch on pressure by pressing «on» – pressure is built up.
  - 2.2 Check flow value.

### Flow decreases:

Thopaz is ready for use.

### Flow does not decrease:

Leak in the system; check tubing connections. Check pressure relief value for tight connection. Replace canister and/ or tubing if necessary.

3. 
  - 3.1 Switch to Standby mode by pressing «Standby» for longer than 3 seconds.

4. Thopaz is in Standby mode and is ready for use with the factory settings.  
(To change factory settings please refer to chapter 11-Changing settings)
5. Connect the Thopaz system to the patients catheter (drain) according to hospital guideline. To start the therapy press «on».

### Safety related tip

The acoustic alarm function «Leak in system – 301» is deactivated for approximately 6 minutes after Thopaz switching on.

## 8.6 Check therapy progress

Two different modes are available. Toggle between the two modes by pressing «Data» or «Graph». In the data mode the set pressure value and the air flow value are digitally displayed. The air flow value is displayed in real time equivalent to the respective parenchymal leakage.

Air flow values are displayed as follows:


- Air flow values between 0 and 1000ml/min are displayed in 10ml steps
- Air flow values greater than 1000 ml/min are displayed in 100ml steps.

### Rounding Function


Digital display of air flow values is rounded to zero ml/min (0ml/min) if the air flow measurement indicates an air flow <5.0ml/min.

If the air flow measurement indicates an air flow  $\geq 5.0$ ml/min and < 15ml/min the digital display of air flow values is rounded to ten ml/min (10ml/min)

## Graph mode


The Graph mode shows the flow and pressure progression as a function of time. The graph is updated every 10 minutes. 4 hours are necessary to have a representative graph. After 60 seconds, the display switches back to the data mode.

The current flow and pressure values are indicated in the information field.


- The left scale shows the pressure in the desired unit (max. -10kPa). The pressure is represented by a filled area.
- The right scale shows the flow in the desired unit (max. 1000ml). The flow values are shown as a line.

## Zoom mode

- 1.1 In Graph mode press [  ] to activate zoom function.

The flow scale is reduced from 1000ml/min. to 100ml/min.


- 2.1 Press [  ] to return to graph function.

If the flow rate does not exceed 100ml/min. for at least 24 hours, Thopaz will display the graph in the zoom mode directly.

## 8.7 Catheter check


- The catheter check is only visible in zoom mode when the flow is 0ml/min.
- **Ask the patient to breathe slowly and deeply 3 times.**
- In zoom mode, a pressure difference from minimal 0.3 kPa between inspiration [  ] and expiration [  ] is displayed at a flow of 0ml/min. This pressure difference is an indication for the respiration cycle of the patient and confirms that the patient catheter (drain) is not clogged.

- Patient catheter (drain) open
- Patient catheter (drain) clogged?


**Important:** The catheter check works only with a patient catheter (drain) in the interpleural space.

## 8.8 Switching Thopaz off

1. 
  - 1.1 Clamp patient tubing with the tube clamp.
  - 1.2 Press «Standby» for longer than 3 seconds to switch off the pressure (Standby mode).
2. 
  - 2.1 Press [  ] quickly to switch off Thopaz.
3. Release, remove and seal the canister by the canister seal. Dispose of canister and patient tubing in accordance with internal hospital guidelines.

## 8.9 Transferring data to PC with ThopEasy+

With ThopEasy+ all stored data can be transferred to a PC as documentation and for completion of the patient records. The data can be expanded with patient information, saved and printed. Operating instructions and the software can be found on the Thopaz CD [REF 079.0035](#) included in the scope of delivery.

For the connection between Thopaz and a PC use the USB cable [REF 079.0034](#).

ThopEasy+ is compatible with the following Microsoft Operating Systems:

Windows 7 32 Bit, Windows 7 64 Bit, Windows XP 32 Bit (SP3 or higher), Windows XP 64 Bit (SP3 or higher), Windows 8 32 Bit, Windows 8 64 Bit, Windows 10.

### WARNINGS

Data transfer via USB is not allowed during therapy.

# 9 Replace canister


## WARNINGS

The canister is replaced on the basis of a visual check or according to the instructions on the display (warning signal). Wear gloves for all operations.

**The warning «canister full» is triggered only once.** After cancelling the warning with the selection buttons, the warning is deactivated and a symbol [  ] flashes in the run mode. After deactivation of the warning, the user is responsible for checking the canister filling level. Please refer to chapter “13 – Troubleshooting and alarms” for a more detailed explanation of the canister full warning.

1. Prepare sterile canister (with/without solidifier).
2. Close patient tubing with the tube clamp.

3.

3.1 Switch to Standby mode by pressing «Standby» for longer than 3 seconds.

4.

4.1 Release the canister by pressing the release button.  
4.2 Remove canister.

5.

5.1 Unpack the canister.  
5.2 Position the openings at the top and attach the bottom part of the canister to Thopaz.


6.

6.1 Push the canister into Thopaz until you hear a click.

**Safety related tip**  
Using a 2l canister refer to chapter 8.

7.  7.1 Press «on». Pressure is built up.

8. Open the tube clamp.
9.  9.1 Check flow value.

10.  10.1 Seal the used canister with the canister seal.

11. Dispose of the used canister in accordance with internal hospital guidelines.

**Canister with solidifier:  
activation of the  
solidification process**

- 
-  a Check if canister is sealed with the canister seal.
- b Press – to open the solidifier chamber.
-  c Shake – to activate the solidification process.
- d Dispose of the used canister in accordance with internal hospital guidelines.

# 10 Changing pressure


## WARNINGS

The settings may only be changed by the physician or on medical orders.

### Change pressure

**Important:** Thopaz is running (Data- or Graph mode).

- 1.1 Press [▶] and [◀] simultaneously. The pressure is indicated on the display.
- 1.2 Set the desired pressure with [▲] and [▼], and confirm with «OK».

### Gravity mode

For patients who are to be treated by gravity drainage (= water seal), the gravity mode can be activated. This mode corresponds to a pressure of  $-0.8\text{ kPa}$ / $-8\text{ cmH}_2\text{O}$ / $-8\text{ mbar}$ / $-6\text{ mmHg}$ . To activate the gravity mode press the [▶] button and confirm with «OK».

# 11 Changing settings


## WARNINGS


The settings may only be changed by the physician or on medical orders.

### Change settings

**Important:** Thopaz must be in Standby mode

1.1 Press [▶] and [◀] simultaneously. The settings are displayed.


- 2.1 Set the desired parameter by pressing [▲] or [▼] and confirm with «OK».
- 2.2 Change the desired parameter by pressing [▲] or [▼] and confirm with «OK».
- 2.3 Press «back» to reactivate the Standby mode.

**Selection**  
 max. -10kPa/-100cmH<sub>2</sub>O  
 kPa/cmH<sub>2</sub>O\*/mbar/mmHg  
 mL/min/L/min  
 off/on  
 no/yes

\* cmH<sub>2</sub>O is not a SI unit according to EU Directive 80/181/EEC.

**The setting is changed for the actual therapy.**


## Factory settings

Thopaz is delivered with the following factory settings. These settings can be changed and saved as the new standard.

Pressure	- 2.0
Unit Pressure	kPa
Unit Flow	mL/min
Leakage	off
Set as standard	no

## Setting a changed parameter as standard

1. Switch to «Set as standard» by pressing [  ] 4 times.
2. Activate «Set as standard» with «OK».
3. Press [  ] or [  ] to change from «no» to «yes» and confirm with «OK». (The display will switch back to «no»).
4. Press «back» to reactivate the Standby mode.


**The new setting is active whenever Thopaz is switched on.**

# 12 Taking a drainage sample


## WARNINGS

To ensure Thopaz can reduce the pressure, it is essential to wait 30 seconds between taking the sample from the sample port and unclamping the patient catheter (drain).


1. Before taking a drainage sample make sure, that there is fluid in the patient tube.
2. Clamp patient catheter (drain).

3.1 Switch to Standby mode by pressing «Standby» for longer than 3 seconds.

\* Syringe 17 G (1,4 mm) or thinner.

4. Disinfect the sample port before taking the sample with a standard disinfectant (e.g. Meliseptol rapid from B.Braun Medical AG).

5.1 Remove with syringe\* air from the patient tubing. Repeat until fluid has gathered at the sample port.

6. With a syringe\* take a sample from the patient tube.

7.1 Switch Thopaz on by pressing «on» – pressure is built up.

8. **Wait 30 seconds!**  
To ensure Thopaz can reduce the pressure, it is essential to wait 30 seconds between taking the sample and unclamping the patient catheter (drain).
9. Unclamp the patient catheter (drain).

# 13 Troubleshooting and alarms

## Pump not running

Check if:

- the mains voltage is correct and if the mains plug is inserted correctly into the fixed mains socket.
- Thopaz is switched on.

## Description

Thopaz distinguishes between warnings, alarms and internal errors. If Thopaz detects any of these situations, an acoustic warning signal sounds and a description of the problem appears on the display. By pressing the two selection buttons simultaneously, the acoustic alarm is suppressed for 60 seconds and the troubleshooting instructions appear on the display.

**Warning**  
Pressure is maintained.

beep

**Alarm**  
Pressure is interrupted.


beep

**Internal error**  
Problem cannot be rectified by the user.  
The error number is displayed.

beep

Visible on display						
Error	Problem description	Troubleshooting instructions	Remarks/Potential source of error	Pressure		
Warning	401	Battery low	Charge battery – connect Thopaz to mains	Remaining time of battery approx. 30min.	Yes	
	306	Canister full	Replace canister	If the canister is not full, see next page	Yes	
	405	Standby mode	Switch Thopaz on or off	In standby mode after 5 minutes	No	
	404	Leakage	<ul style="list-style-type: none"> <li>– Check patient status</li> <li>– Check system for leaks</li> <li>– Continue with 'OK'</li> </ul>	<ul style="list-style-type: none"> <li>– Irregularity in the course of therapy</li> <li>– Indication of leakage</li> <li>– See Chapter 6</li> </ul>	Yes	
	402	USB connection not permitted	Unplug USB cable	USB connection not permitted during operation or in the vicinity of the patient	Yes	
	406	Temperature too high	Prepare backup device	Do not expose Thopaz to heat sources	Yes	
Alarm	305	Battery empty	Charge battery – connect Thopaz to mains		No	
	302 304 308	System clogged	<ul style="list-style-type: none"> <li>– Check that tubing is clear/not twisted</li> <li>– Replace canister</li> </ul>	<ul style="list-style-type: none"> <li>– Tubing is kinked or clogged, or filter in canister is clogged for approx. 10 min.</li> <li>– Thopaz detects clotting of the Thopaz tubing, NOT of the patient catheter</li> </ul>	No	
	301	Leak in system	Check system for leaks	<ul style="list-style-type: none"> <li>– Disconnection</li> <li>– Missing seal between Thopaz and canister</li> </ul>	No	
	311 314	Selftest failed	Snap the canister out and in again	Do not switch Thopaz on if the patient is already connected	No	
	313	Filter clogged	<ul style="list-style-type: none"> <li>– Replace canister</li> <li>– Continue with 'on'</li> </ul>	Filter in canister clogged	No	
	315	Thopaz overheated	Replace Thopaz		No	
	Internal error		Internal error XXX	<ol style="list-style-type: none"> <li>1. Switch Thopaz off and on again</li> <li>2. Inform Medela Customer Service</li> </ol>	If error message repeats, inform Medela Customer Service	

### Safety related tip

Warning «canister full» / Alarm «filter clogged»

The warning «canister full» is triggered when the canister filling level approximately reaches the maximal scale line indicated on the respective canister label. Change the canister according the instructions for use, chapter 9 – Replace canister. If secretions line the walls of the canister, the warning «canister full» may be activated prematurely.

The warning «canister full» is triggered only once. After cancelling the warning with the selection buttons, it is deactivated and a symbol [  ] flashes in the run mode. After deactivation of the warning, the user is responsible for checking the canister filling level.

The warning is reactivated under the following conditions:

- if Thopaz is switched to standby and then on again – for example, on replacing the canister.
- if the canister full sensor is not activated for at least 10 minutes.

An overfilling of the respective canister volume will trigger the alarm «filter clogged». The alarm is triggered due to the clogging of the hydrophilic filter upon contact with liquids or secretions. Please refer to chapter «5 – Overview» for location of the hydrophilic, overflow protection/bacteria filter in the canisters. Subsequent to a clogged filter and the alarm «filter clogged» the vacuum applied to the patient is interrupted and cannot be regulated or maintained. An additional accumulation of liquids, secretions and air is released through the pressure relieve valve at positive pressures of 0.2 – 0.5 kPa. This pressure release function is not a substitute for a free overflow function. Therefore the canister must be replaced once full. Please refer to chapter «5 – Overview» for location of the pressure relieve valve. Change the canister according the instructions of use, chapter «9 – Replace canister».

### Safety related tip

The acoustic alarm «Leak in system – 301» is deactivated for approximately 6 minutes after switching Thopaz on.

# 14 Warranty, maintenance and checks

## **Warranty**

Medela AG warrants the device will be free from defects in materials and workmanship for a period of 2 years from the date of delivery ex works. Faulty material will be replaced free of charge during this period if not resulting from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with our appliances.

In no event shall Medela AG be liable for claims which exceed the scope of warranty described including liability for consequential damages, etc. The right to the replacement of faulty parts will not be recognized by Medela if any work has been carried out on Thopaz by unauthorized persons. This warranty is subject to the appliance being returned to a Medela service centre.

## **Maintenance**

Thopaz is a maintenance free device. Each time the device is switched on, a self-test is carried out which checks the internal functions of Thopaz. A «beep» indicates that Thopaz successfully passed the self-test and that the device is switched on. Safe operation is ensured by the functional check including canister and tubing to be carried out before each use.

The safety instructions must be observed.

## **Routine Check**

Medela recommends subjecting Thopaz to a routine check once a year. The work to be carried out and a protocol template in accordance to EN/IEC 62353:2008 is given in Appendix B in these instructions.

## Instructions for FPT test


**Important:** Patient must not be connected! Thopaz must be disconnected from the mains (battery operation) and turned off.

(Tube e.g.  $\varnothing 12 \times 7$  mm, not included in the scope of delivery) 0.8l

1. Connect the tubing (079.0021) and the 0.8l (079.0016) canister to Thopaz according to chapter 8.
2. Connect vacuum gauge with tube to the patient connector.

3.

Press simultaneously

- 3.1 Press and hold [▶], [◀] simultaneously, then...
- 3.2 Press additionally [⏻] to start safety test.
- 3.3 Follow the instructions on the display.

## Safety Check

There are no prescribed or recommended safety checks to be carried out.

### Rationales

- The engineering of the Medela Thopaz has been verified by independent test institutes for compliance with the EN / IEC 60601-1 standard. Copies of such certificates may be obtained on request. Because of the method of construction used, Medela does not expect electrical safety to be affected at any time during the life of the product - provided that Thopaz is repaired throughout its service life strictly and exclusively by Medela authorized service centers and that it is used properly in accordance with the intended use.
- For Thopaz as a device in protection class II (EN IEC 60601-1), there is no protective earth conductor; protective earth related tests (protective earth resistance, earth leakage current etc.) are not applicable.
- Thopaz enclosure is made entirely of insulating material. Tests of the touch current using common measuring instruments will therefore not reveal measurable values.
- The applied parts of Thopaz are connected to Thopaz via non-conductive vacuum hoses, canisters and tubings. Even when suctioning a conductive fluid (9 g/L NaCl; EN ISO 10079-1;1999) until the overflow protection device activates, measurements of the patient leakage current using common measuring instruments will not reveal measurable values.
- Thopaz systems do not have patient connections and functional earth conductors; tests related to those connections and currents are not applicable.

# 15 General cleaning guidelines

## **WARNINGS**

After each use, the parts that have been in contact with aspirated secretions are to be cleaned and disinfected or disposed.

## **General notes**

- For detailed information refer to Cleaning & Reprocessing Guide: 200.2391.
- This cleaning guideline relates to the cleaning of the entire Thopaz product ranges see Appendix A.
- Never place electrical devices in water or other liquids, specifically Thopaz and the docking station.
- Do not spray or pour liquid directly onto Thopaz/docking station.
- Thopaz product range cannot be sterilized.
- Rinsing the Thopaz product range in a washing machine is not permitted.
- Immersion disinfection, thermal disinfection and ultrasound cleaning for the Thopaz product range are not permitted.
- Higher temperature than 45°C, 113°F can cause protein coagulation which can lead to problems further in the process.
- Clean all surfaces immediately after use to avoid residues from drying and to prevent growth of microorganisms.
- Do not use cleaning agent/detergent based on phenol.
- Do not use steel brushes or steel wool for cleaning.
- Store medical products dry and dust free.

## **Disposables**

This symbol indicates a single use product. This product is not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

---

## 1. Thorough cleaning of Thopaz

Thorough cleaning can either be done at the point of use with sterile cold water only or in dedicated health care facilities with the additional use of enzymatic detergents according the manufacturer's instructions for use.

- 1.1 Disconnect the power plug from the power socket.
- 1.2 Disinfect your hands and put on disposable gloves and proper protective gear.
- 1.3 Separate all individual parts. Dispose single use products in accordance with internal hospital guidelines.
- 1.4 Remove visible dirt with lint free nonwoven wipes wetted with sterile cold water (<40 °C, 104 °F) Clean around the vacuum port to ensure neither sterile water nor visible dirt enters the vacuum port.
- 1.5 Ensure that all surfaces are thoroughly wetted and keep moist for a minimum of 5 minutes then repeat this step with another wipe.
- 1.6 Dispose protective gloves and disinfect your hands.

---

## 2. Intermediate level disinfection

- 2.1 Disinfect your hands and put on disposable gloves and proper protective gear.
- 2.2 Use disinfecting wipes\* according the manufacturer's instructions for use. Clean around the vacuum port to ensure neither disinfecting agent nor visible dirt enters the vacuum port.
- 2.3 Let product work in for 5 minutes and wipe afterwards with another disinfecting wipe\*.
- 2.4 Allow the surface to dry for a minimum of 5 minutes.
- 2.5 Dispose protective gloves and disinfect your hands.

\*Recommended agent for intermediate level disinfection:

CaviWipes®

Metrex® Research

Address: 1717 West Collins Avenue

Orange, CA 92867, U.S.A.

Homepage: <http://www.metrex.com>

Phone: (800) 841 1428

Email: [metrexcustcare@sybrondental.com](mailto:metrexcustcare@sybrondental.com)

Mikrozid® AF Wipes

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Address: Robert-Koch Str. 2,

22851 Norderstedt, GERMANY

Homepage: <http://www.schuelke.com>

Phone: +49 (0) 40 521 00 0

Email: [info@schuelke.com](mailto:info@schuelke.com)

---

## 3. Storage

- 3.1 After appropriate cleaning and disinfecting check the device before re-use.
  - 3.2 Store products dry and dust free under storing conditions (see chapter 18 – Technical specifications).
-

# 16 Accessories

## **WARNINGS**

Thopaz was verified in combination with the accessories listed in Appendix A. For a correct and safe operation use Thopaz with these accessories only. Further information is supplied on the instruction sheet of the individual accessory. For Accessories overview/legend refer to Appendix A.

# 17 Signs and symbols

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates a class II device.

This symbol indicates the date of manufacture.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates a type CF applied part.

This symbol indicates that the device should not be used after the date shown.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates the legal specifications of the system.

This symbol indicates a prescription device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. (for US only)

This symbol indicates manufacturer's catalogue number.

This symbol indicates to follow instructions for use.

This symbol indicates MR UNSAFE

This symbol indicates manufacturer's serial number.

This symbol indicates manufacturer's batch code.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates that interferences may occur in the vicinity of equipment marked with this symbol.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates safety related tip.

This symbol indicates the manufacturer.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the class of the system.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates to handle the fragile device with care.

This symbol indicates that device contains phthalates. Specific phthalates contained in the device are identified below the symbol. DEHP: Di(2-ethylhexyl) phthalate

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

This symbol indicates the maximum vacuum level of the system.

This symbol indicates that battery is empty.

This symbol indicates do not use the device if package is damaged.

This symbol indicates that battery is almost empty.

This symbol indicates the flow levels of the system.

This symbol indicates that battery is being charged (moving bars).

## unit

This symbol indicates a specified amount of quantity.

This symbol indicates the electrical specifications of the system.

This symbol indicates that battery is fully charged.

This symbol indicates to not dispose of the device together with unsorted municipal waste (for EU only).

This symbol indicates alternating current.

This symbol indicates to turn system on into standby mode.

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates direct current.

This symbol indicates a notification.

This symbol indicates a carton package.

This symbol indicates the weight of the system.

This symbol indicates that audio is paused.

This symbol indicates to keep the device away from sunlight.

This symbol indicates the dimensions (h x w x d) of the system.

This symbol indicates up.

## pcs

This symbol indicates number of items.

This symbol indicates down.

# 18 Technical specifications

low vacuum,  
- 10 kPa / - 75 mmHg /  
100 cmH<sub>2</sub>O/mbar  
Tolerance: +/- 15 %

Measured at 0m, atmospheric pressure: 1013.25 hPa Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

5 l/min. (+/- 10%)

1,0 kg / 2.2 lbs

Power supply  
Model: TR30RAM120  
IEC: 60601-1  
Input: 100-240V~,  
0.8-0.4 A,  
47-63 Hz  
Output: 12V, 2.5A

Pump  
12 VDC, 20W

ISO 13485  
CE (93/42/EEC), IIb


0123

H x W x D  
223 x 255 x 95 mm

Transport/Storage Conditions

Operating Conditions

**IP33**

# 19 Disposal

Thopaz comprises metals and plastics and should be disposed of in accordance with the European directives 2011/65/EC and 2012/19/EC. Additional, local guidelines must also be observed. The electronic components must be disposed of separately, in accordance with the local regulations. This product contains lithium-ion batteries which bear risk of fire, explosion and burns, if disposed of improperly. Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire. Dispose of in accordance with your local requirements.

Please take care that you dispose of Thopaz and its accessories in accordance with the hospital disposal guidelines.

## **User information for the disposal of electrical and electronic equipment**

This symbol means that the electrical and electronic equipment must not be disposed as normal household refuse. A correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about the disposal contact the manufacturer, your local caregiver or healthcare provider. This symbol is only valid in the European Union. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.

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# 20 Technical documentation

## **WARNINGS**

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of Thopaz device.

## **EMC**

Thopaz is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfill the requirements of the relevant IEC 60601-1-2 standard. Untested HF (high-frequency) sources, radio networks or the like can influence the operation of the device and may not be operated in combination with the system. Thopaz is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information. Portable and mobile RF communication devices (mobile telephones) can affect Thopaz.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007, Table 1)

### Electromagnetic emissions

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that they are used in such environment.

Emission Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thopaz uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Thopaz is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

#### WARNINGS

Thopaz should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Thopaz should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007, Table 2)

### Electromagnetic immunity

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 4 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Thopaz requires continued operation during power mains interruptions, it is recommended that Thopaz be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

**NOTE**  $U_T$  is the a.c. mains voltage prior to application of the test level.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007, Table 4)

### Electromagnetic immunity

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that Thopaz is used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of Thopaz, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 0.35 \sqrt{P}$ $d = 0.35 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 0.7 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Thopaz are used exceeds the applicable RF compliance level above, Thopaz should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Thopaz.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

### Recommended separation distance between portable and mobile RF communications equipment and Thopaz

Thopaz is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Thopaz can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Thopaz as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter	Separation distance according to frequency of transmitter m (meters)		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
W			
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.21
100	3.50	3.50	7.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# Appendix A

**REF**

	079.0000	Thopaz – medical aspirator for thoracic drainage with Instructions for use 200.0685
	079.0002	Thopaz – medical aspirator for thoracic drainage with Instructions for use 200.0686
	079.0003	Thopaz – medical aspirator for thoracic drainage with Instructions for use 200.2004 / 101034234
<b>Disposables</b>	079.0011	Thopaz canister 0.3 L, antifoaming*
	079.0012	Thopaz canister 0.3 L with solidifier, antifoaming*
	079.0016	Thopaz canister 0.8 L, antifoaming*
	079.0017	Thopaz canister 0.8 L with solidifier, antifoaming*
	079.0018	Thopaz canister 2 L, antifoaming*
	079.0019	Thopaz canister 2 L with solidifier, antifoaming*
	079.0021	Thopaz tubing, single*
	079.0022	Thopaz tubing, double*
	079.0023	Thopaz tubing, single, small connector*
	079.0024	Thopaz tubing, double, small connector*
	079.0025	Thopaz tubing, single, large connector*
	079.0026	Thopaz tubing, double, large connector*
<b>Accessories</b>	079.0031	Carrying strap
	079.0037	Docking station
	079.0038	Adapter docking station for Thopaz canister 2 L
	079.0034	USB cable
	079.0036	Holder with standard rail
	079.0039	Thopaz sealing cap*
	079.0040	Universal holder with standard rail
	077.1456	Vacuum gauge
<b>Documents</b>	079.0035	Thopaz CD
	200.0685	Thopaz Instr. for use (EN/DE/FR/IT/NL/SV/DA/ES/PT/PL)
	200.0686	Thopaz Instr. for use (FI/NO/CZ/SK/HU/EL/LV/LT/ET/TR)
	200.2004	Thopaz Instr. for use (SL/HR/BG/RO/SR/RU/KO/JP/AR)
	101034234	Thopaz Instr. for use (EN/ZH)
<b>Spare parts</b>	077.0148	Thopaz wall mains adapter coaxial DC plug
	077.0117	Screw covers
	077.0119	Feet
	099.0343	Seal orange
	077.0124	Cover ports
	077.0129	Set of holders with tool
	099.0348	Release button

\* **STERILE EO**

# Appendix B

DIN EN 62353: 2008-8

<b>Hospital</b>	Thopaz serial number:
	Firmware version:
Location:	Run time (Settings 2/2): h

Visual checks	NO	YES
Cracks, brittle areas, discolouration on housing and display		
Covers for ports are present and can be clicked in		
Standard rail and bed holders are present and can be opened		
Seal is present, is clean and not brittle		
Release button can be pressed, canister can be clicked in		
Thopaz is clean with no dirty areas		

Tests	OK	Not OK
<b>Safety test</b> Carry out test as described in Instructions for use. Thopaz must be disconnected from the mains (battery operation). Required test equipment: vacuum gauge -0-10kPa Class 2.5 % or better (E.G.077.1456)	<input type="checkbox"/>	<input type="checkbox"/>
<b>Acoustic signal (beep)</b> Audible?	<input type="checkbox"/>	<input type="checkbox"/>
<b>Battery</b> Connect Thopaz to mains adapter or place it in docking station.	<input type="checkbox"/> Battery is being charged.	<input type="checkbox"/> Battery is not being charged.

Overall evaluation	Errors	Actions
No safety or functional technical errors were detected.	<input type="checkbox"/>	Thopaz is ready for use.
Errors that are not safety relevant were detected. No direct risk.	<input type="checkbox"/> 1 _____ 2 _____ 3 _____	Thopaz is ready for use, errors must be corrected by __/__/__
Safety relevant errors were detected.	<input type="checkbox"/> 1 _____ 2 _____ 3 _____	Do not use Thopaz until the errors are corrected. Send Thopaz to Medela Service Centre.

Next test: \_\_/\_\_/\_\_\_\_\_

Tester/Signature:	Place/Date:
-------------------	-------------

# Notes

# Notes

Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland  
www.medela.com

**International Sales**

Medela AG  
Lättichstrasse 4b  
6341 Baar  
Switzerland  
Phone +41 41 562 51 51  
Fax +41 41 562 51 00  
ism@medela.ch  
www.medela.com

**Australia**

Medela Australia Pty Ltd,  
Medical Technology  
3 Arco Lane  
Heatherton, Vic 3202  
Australia  
Phone +61 3 9552 8600  
Fax +61 3 9552 8699  
contact@medela.com.au  
www.medela.com.au

**Canada**

Medela Canada Inc.  
4160 Sladeview Crescent Unit # 8  
Mississauga, Ontario  
L5L 0A1  
Canada  
Phone +1 905 608 7272  
Fax +1 905 608 8720  
info@medela.ca  
www.medela.ca

**France**

Medela France SAS  
14, rue de la Butte Cordière  
91154 Etampes cedex  
France  
Phone +33 1 69 16 10 30  
Fax +33 1 69 16 10 32  
info@medela.fr  
www.medela.fr

**Germany**

Medela Medizintechnik  
GmbH & Co. Handels KG  
Postfach 1148  
85378 Eching  
Germany  
Phone +49 89 31 97 59-0  
Fax +49 89 31 97 59 99  
info@medela.de  
www.medela.de

**India**

Medela India Private Limited  
602-603 | 6th Floor,  
JMD Megapolis | Sohna Road,  
Sector – 48 | Gurgaon 122 001  
India  
Phone +91 124 4416999  
Fax +91 124 4416990  
info@medela.in  
www.medela.in

**Italy**

Medela Italia Srl  
Via Turrini, 13/15 – Loc. Bargellino  
40012 Calderara di Reno (BO)  
Italy  
Phone +39 051 72 76 88  
Fax +39 051 72 76 89  
info@medela.it  
www.medela.it

**Japan**

Medela K.K.  
G-PLACE Sasazuka 13F  
2-26-2, Sasazuka, Shibuya-ku  
Tokyo 151-0073  
Japan  
Phone +81 3 3373 3454  
Fax +81 3 3373 3457  
info@medela.jp  
www.medela.jp

**Netherlands & Belgium**

Medela Benelux BV  
Uilenwaard 31  
5236 WB's-Hertogenbosch  
Netherlands  
Phone +31 73 690 40 40  
Fax +31 73 690 40 44  
info@medela.nl  
info@medela.be  
www.medela.nl  
www.medela.be

**Poland**

Medela Polska Sp. z o.o.  
Ul. Lewinowska 8  
03-684 Warszawa  
Poland  
Phone +48 22 865 12 50  
Fax +48 22 864 38 65  
info@medela.pl  
www.medela.pl

**Russia**

OOO Medela  
Vavilova str. 97, entrance 3  
117335 Moscow  
Russia  
Phone +7 495 980 6194  
Fax +7 495 980 6194  
info@medela-russia.ru  
www.medela-russia.ru

**Spain & Portugal**

Productos Medicinales Medela, S.L.  
C/ Llacuna, 22 4ª Planta  
08005, Barcelona  
Spain  
Phone +34 93 320 59 69  
Fax +34 93 320 55 31  
info@medela.es  
info@medela.pt  
www.medela.es  
www.medela.pt

**Sweden**

Medela Medical AB  
Box 7266  
187 14 Täby  
Sweden  
Phone +46 8 588 03 200  
Fax +46 8 588 03 299  
info@medela.se  
www.medela.se

**Switzerland**

Medela AG  
Lättichstrasse 7  
6341 Baar  
Switzerland  
Phone +41 848 633 352  
Fax +41 41 562 14 42  
contact@medela.ch  
www.medela.ch

**United Kingdom**

Medela UK Ltd.  
Huntsman Drive  
Northbank Industrial Park  
Irlam, Manchester M44 5EG  
UK  
Phone +44 161 776 0400  
Fax +44 161 776 0444  
info@medela.co.uk  
www.medela.co.uk

**USA**

Medela LLC  
1101 Corporate Drive  
McHenry, IL 60050  
USA  
Phone +1 877 735 1626  
Fax +1 815 307 8942  
info-healthcare@medela.com  
www.medela-healthcare.us