CAESAID CUP

CAESAID CUP

This obstetrical instrument is utilized to assist a clinician in the delivery of an infant during caesarean section.

Important note

Proper surgical procedures and techniques are the responsibility of the medical practitioner. Vacuum assisted delivery must only be performed or supervised by trained and experienced clinicians. Each physician must evaluate the suitability of the procedure based on his or her own medical training and experience.

General information

The instrument consists of a 60mm vacuum cup with a moulded collar and integral vacuum tube. Manufactured in latex-free silicone, it is flexible and easily compressed for insertion. The instrument should be connected to an electric vacuum pump of the type used in vaginal vacuum delivery and fitted with an external reservoir jar.

Pump requirements

- It is not recommended for use with a hand held vacuum pump, since liquor may be aspirated during establishing a vacuum
- This obstetrical instrument has been developed for use in conjunction with an external vacuum source capable of precise vacuum regulation to ensure proper vacuum levels (-60 to -80kPa/-4.50 to -6.00mmHg)
- The vacuum source should be an electric pump designed for medical purposes.
- We recommend using the Medela Basic suction pump.
- Note: Wall suction is not to be used.

Tubing requirements

The tubing connecting the instrument to the vacuum source is vital for safe and proper function of the instrument. To ensure this only Medela original tubing should be used.

Cleaning/Disinfection /before use

After use, soak the suction cups and tubing at room temperature for 30 minutes in a disinfectant (cleaning solution A).

- After the reaction time, rinse the residue of the solution from the parts with water. Rinse the tubing.
- Clean suction cup and tubings separately with cleaning solution B.
- Use a brush to remove residue.
- Rinse the tubing for at least 2 minutes with the cleaning solution.
- The parts should be free from residue and foreign particles.
- Rinse the parts and the tubing with water (preferable fully desalinated) for 2 minutes.
- Then dry the parts with a fibrous cloth.
- Hang the tubing overnight to dry.

A) Fluid aldehyde-free disinfection solution with cleaning effect for use as a soaking bath, bactericidal, fungicidal, limited virucidal, with good material compatibility for stainless steel, non-ferrous material and plastics including silicone. For the validation of the procedure a 1% neodesher®MediMym-solution is used. For the validation of the procedure a 2% neodesher®Septo MED-solution is used.

Caution

Although the use of vacuum-aided-delivery (VAD) in vaginal delivery is women living with HIV has been studied and found to be relatively safe, the use of the CaesAid VAD to caesarean delivery in women living with HIV has not been assessed. For this reason, CaesAid VAD in HIV positive women at caesarean section cannot be recommended at this stage.

Disposal

Dispose of according to local regulations

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

The use of other cleaning and disinfecting agents should be cleared with the hospital hygiene administrator.

Sterilization (also before first use)

Suction cups and tubings must be assembled and packed before sterilization.

Sterilization: In the autoclave with saturated steam at 134 °C, with triple pre-evacuation, sterilization time is 6 minutes.

CE 0123

Made in Switzerland

Medela AG
Lättichstrasse 4b,
6340 Baar, Switzerland
www.medelahealthcare.com

International Sales
Medela AG
Lättichstrasse 4b, 6340 Baar
Phone +41 41 562 51 51, Fax +41 41 562 51 00
customercare@medela.ch, www.medelahealthcare.com

Germany
Medela Medizintechnik GmbH & Co. Handels KG, Postfach 1148, 85378 Eching
Phone +49 89 31 97 59-0, Fax +49 89 31 97 59 99
info@medela.de, www.medelahealthcare.de

Switzerland
Medela AG
Lättichstrasse 4b, 6340 Baar, Switzerland
Phone +41 848 633 352, Fax +41 41 562 51 00
customercare@medela.ch, www.medelahealthcare.ch

ENGLISH

CAESAID CUP

This obstetrical instrument is utilized to assist a clinician in the delivery of an infant during caesarean section.

Expected product life
60 cleaning/sterilisation cycles

Safety rules

- Do not use it:
 - If the device is damaged.
 - If the device shows obvious safety defects.
 - If damaged – especially tears – are noted, the parts may no longer be used.
 - If the CaesAid Cup is discoloured it should be discarded since this may indicate material deterioration.
- Wall suction is never to be used.
- Vacuum extraction must only be performed or supervised by trained and experienced clinicians.
- The neonatal staff must be aware of the indications, contraindications and precautions for assisted delivery and must watch for the signs of complications of vacuum assisted delivery. The CaesAid Cup and VAD system is not recommended for use with classical uterine incisions.
- Ensure against intrusion of extraneous tissues and ensure proper placement of cup.
- Attention: Never apply the cup to any portion of the infant's face as this can cause fetal trauma.
- Insertion should be performed carefully to avoid malfunction, maternal or fetal trauma.
- Discontinue vacuum extraction if the vacuum cup becomes disengaged or pops off twice.
- Parents should be instructed to check the baby's health.
- The infant has to be watched for several days for possible problems.

General information

The instrument consists of a 60mm vacuum cup with a moulded collar and integral vacuum tube. Manufactured in latex-free silicone, it is flexible and easily compressed for insertion. The instrument should be connected to an electric vacuum pump of the type used in vaginal vacuum delivery and fitted with an external reservoir jar.

Pump requirements

- It is not recommended for use with a hand held vacuum pump, since liquor may be aspirated during establishing a vacuum
- This obstetrical instrument has been developed for use in conjunction with an external vacuum source capable of precise vacuum regulation to ensure proper vacuum levels (-60 to -80kPa/-4.50 to -6.00mmHg)
- The vacuum source should be an electric pump designed for medical purposes.
- We recommend using the Medela Basic suction pump.
- Note: Wall suction is not to be used.

Tubing requirements

The tubing connecting the instrument to the vacuum source is vital for safe and proper function of the instrument. To ensure this only Medela original tubing should be used.

Cleaning/Disinfection /before use

After use, soak the suction cups and tubing at room temperature for 30 minutes in a disinfectant (cleaning solution A).

- After the reaction time, rinse the residue of the solution from the parts with water. Rinse the tubing.
- Clean suction cup and tubings separately with cleaning solution B.
- Use a brush to remove residue.
- Rinse the tubing for at least 2 minutes with the cleaning solution.
- The parts should be free from residue and foreign particles.
- Rinse the parts and the tubing with water (preferable fully desalinated) for 2 minutes.
- Then dry the parts with a fibrous cloth.
- Hang the tubing overnight to dry.

A) Fluid aldehyde-free disinfection solution with cleaning effect for use as a soaking bath, bactericidal, fungicidal, limited virucidal, with good material compatibility for stainless steel, non-ferrous material and plastics including silicone. For the validation of the procedure a 1% neodesher®MediMym-solution is used. For the validation of the procedure a 2% neodesher®Septo MED-solution is used.

Caution

Although the use of vacuum-aided-delivery (VAD) in vaginal delivery is women living with HIV has been studied and found to be relatively safe, the use of the CaesAid VAD to caesarean delivery in women living with HIV has not been assessed. For this reason, CaesAid VAD in HIV positive women at caesarean section cannot be recommended at this stage.

Disposal

Dispose of according to local regulations

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

The use of other cleaning and disinfecting agents should be cleared with the hospital hygiene administrator.

Sterilization (also before first use)

Suction cups and tubings must be assembled and packed before sterilization.

Sterilization: In the autoclave with saturated steam at 134 °C, with triple pre-evacuation, sterilization time is 6 minutes.

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Summary- Use of the Caesaid at Caesarean Section

Dispose of according to local regulations

The device may be used for both elective and emergency caesarean section.

Use of the CaesAid cup allows quick, atraumatic delivery of the fetal head.

It possibly allows a smaller incision with less blood loss, and consequently quicker repair.

The uterus is opened in the conventional manner, and the fetal presenting part displayed.

The cup is applied to the posterior fontanelle after excess liquor has been suctioned.

The finger explores the rim of the cup to ensure no tissue or cord is trapped between the cup and the fetal scalp.

The vacuum is immediately increased to –60 kPa by means of an electric pump.

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