

# Medela Vario18 suction pump

Pump Specifications			
Item numbers	<b>AC/DC version</b>		
	026.2110	230–240V / 50Hz	Euro plug 2xA, 1xB
	026.2112	230–240V / 50Hz	Euro plug A, B
	026.2140	2230–240V / 50Hz	UK plug A
	026.2150	230–240V / 50Hz	AUS/NZ plug A
	026.2120	100–120V / 60Hz	US plug A
	026.2160	230–240V / 50Hz	Euro plug A
	<b>AC/DC version</b>		
	026.5110	100–240V / 50–60Hz	Euro plug 2xA, 1xB
	026.5112	100–240V / 50–60Hz	Euro plug A,B
	026.5140	100–240V / 50–60Hz	UK plug A
	026.5150	100–240V / 50–60Hz	AUS/NZ plug A
	026.5120	100–240V / 50–60Hz	US plug A
	Instructions for use	177.2604	EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (A)
	177.2605	ET, HU, CS, LT, LV, PL, SL, SK, DA, EL (B)	
Service manual	177.2610	EN, DE, FR	
Max. vacuum [kPa/mmHg]	-75/-563	Tolerance: ± 15%	
Flow [l/min]	18l/min.	± 10%	
Dimensions [mm]	380 x 170 x 285 height including filter and 90° coupling piece: 440		
Weight [kg]	AC 3.5	AC/DC 4.2	
Drive unit	QuatroFlex piston/cylinder technology		
Power supply	AC	230–240V, 50Hz / 230–240V, 60Hz / 120V, 60Hz	
	AC/DC	100–240V, 50–60Hz	
Power consumption [W]	AC	90 (70 for 100–120V version AC)	
	AC/DC	80	
Noise level [dB(A)]; measured in Medela acoustic room at 1 m distance	51.2		
Vacuum regulator	Membrane		
Protection against electric shock	Type CF (cardiac floating: max. earth leaking current 10 µA)		
Medical Device Classification – According to 93/42/EEC	Class IIa		
Protection against electric shock	Class II Equipment		
Ingress Protection	IP21		
Battery	NiMH		
Operation from battery (at full charge)	30 minutes		
Charging time of battery	5 hours		
Low battery indication	LED and acoustic		
Overflow protection	Disposable overflow protection/bacterial filter (5 pcs included)		
Warranty [years]	2; QuatroFlex and battery 6 months		
Operating temperature, relative humidity, pressure	5 to 40 °C, 15–93 %, 70–106kPa		
Transport/storage temperature, relative humidity, pressure	-20 to +50 °C, 15–93 %, 70–106kPa		
Made in	Switzerland		

## Main compatible accessories

Reusable Collection System	Item number
Reusable 1 l set with PSU suction jar	077.0715
Reusable 2 l set with PSU suction jar	077.0716

Disposable Collection System	Item number
Adapter for Vario	077.0105
Suction jar PC, 1.5L	077.0082
Disposable suction liner 1.5L	077.0083
Disposable suction liner 1.5L with solidifier	077.0084
Suction jar PC, 2.5L	077.0085
Disposable suction liner 2.5L	077.0086
Disposable suction liner 2.5L with solidifier	077.0087

Filters	Item number
Vario disposable filter set	077.0834
Disposable overflow protection/bacteria filter with Medela connections	077.0571
Disposable overflow protection/bacteria filter with Medela and conical connections	077.0577

Tubing, silicone	Item number
Silicone tubing Ø 7x12mm, without coupling pieces 1 m	077.0901
Silicone tubing Ø 7x12mm, without coupling pieces, 2 m	077.0902
Silicone tubing Ø 7x12mm with 1 coupling piece, 60 cm	077.0912
Silicone tubing Ø 7x12mm with 1 coupling piece, 100 cm	077.0913
Silicone tubing Ø 7x12mm with 2 coupling pieces, 30 cm	077.0921
Silicone tubing Ø 7x12mm with 2 coupling pieces, 60 cm	077.0922
Silicone tubing Ø 7x12mm, 25 cm, with 1 coupling piece and 1 coupling piece 90°	077.0931

Tubing, disposable	Item number
Disposable (PVC) tubing, 180 cm, with fingertip, sterile	077.0951
PVC tubing Ø 12.7x18mm, 200 cm for patient connection Ø 10–14 mm	077.0963
Disposable PVC tubing, 150 cm CH 25, sterile	077.0965

Optional	Item number
Clamp holder, plastic	077.0521
Wall bracket, single	077.0523
Trolley small	077.0511
Vario rail holder 10x25 mm	077.0821
Vario universal holder	077.0823
Vario carrying bag	077.0831
Car connection cable for Vario 12V DC	077.0832
Disposable specimen cup, Ø 6–10 mm (for reusable collection system)	077.0562
Disposable specimen cup for disposable collection system	077.0094

CAUTION: Federal US law restricts this device to sale by or on the order of a physician.

Medical Vacuum Technology  
for Healthcare Professionals

Please contact us or your local Medela  
representative for details.

 Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland  
[www.medela.com](http://www.medela.com)

**Local contact:**